

TINGLE X-RAY PRODUCTS, Inc.

TECHNIQUE CHART HIGH FREQUENCY

200 SPEED

	POS.	KVP	DIST	GRID	FOCAL	CORRECTION																						
THICKNESS OF PART IN CENTIMETERS						2	4	6	8	10	12	14	16	18	20	22	24	26	28	30	32	34	36	38	40	42	44	46
CERVICAL SPINE	A.P.	70	40	YES	LARGE			4	6	12	14	20	24	36	40	60	70	100	120	160								
CERVICAL SPINE	LAT.	80	72	YES	LARGE			10	16	20	24	32	40	60	80	120	150	180	240	300								
THORACIC SPINE	A.P.	80	40	YES	LARGE					10	12	20	20	30	36	50	60	70	100	120	150	160	200					
THORACIC SPINE	LAT.	80	40	YES	LARGE					6	10	12	20	20	24	30	36	40	60	70	80	120	150	200	300	360		
LUMBAR SPINE	A.P.	90	40	YES	LARGE							6	10	16	24	30	40	60	80	120	140	180	240	300	360	400		
LUMBAR SPINE	LAT.	90	40	YES	LARGE							20	24	24	30	40	50	60	70	80	120	140	160	200	300	360		
PELVIS&HIP	A.P.,FROGLEG	80	40	YES	LARGE					6	12	20	24	36	40	70	80	120	160	220	360	400						
PELVIS&HIP	TRUE LAT	80	40	YES	LARGE					10	14	20	20	24	30	36	40	60	70	80	80	100	120	160	200			
PELVIS&HIP	PELVIMETRY	90	40	YES	LARGE							24	24	36	40	50	60	70	80	120	140	160	200	240				
THICKNESS OF PART IN CENTIMETER						2	4	6	8	10	12	14	16	18	20	22	24	26	28	30	32	34	36	38	40	42	44	46
ABDOMEN	P.A.	80	40	YES	LARGE					6	10	16	24	36	50	70	100	160	180	240	360	400						
RIBS	P.A. AD	70	40	YES	LARGE					6	10	12	14	20	20	24	30	36	40	50	60	70	80	100	120			
RIBS	P.A. BD	85	40	YES	LARGE					10	14	16	20	24	30	40	50	70	40	100	120	160	200	220	300			
CHEST	P.A.	110	72	YES	LARGE									2	3	3	4	4	6	6	8	10	16	20	24	30		
CHEST	LAT.	110	72	YES	LARGE									2	3	4	4	6	6	8	10	12	16	20	24	36	50	70
THICKNESS OF PART IN CENTIMETERS						2	4	6	8	10	12	14	16	18	20	22	24	26	28	30	32	34	36	38	40	42	44	46
EXTREMITIES	ALL	40	40	NO	SMALL			3	6	8	12	20	24	30	36	50												
EXTREMITIES	ALL 65	40	YES	SMALL				2	4	6	10	12	14	20	24	36	40	50	60	70	80							
SKULL	LAT	75	40	YES	LARGE						10	14	20	24	30	36	30	40	50	70	80							
MANDIBLE, SKULL	P.A.	75	4-	YES	LARGE				6	10	14	24	36	40	50	60	70	80	100									
MANDIBLE	LAT., OBL	40	40	NO	SMALL			3	6	8	12	20	24	30	36	50												
SINUS(NASAL)	LAT	70	40	YES	LARGE						14	20	24	36	40	50	60	70	80	100								
WATER,CALDWELL,TMS.LAW		75	40	YES	LARGE					14	20	24	36	40	60	70	80	100	120	160								